


OUR HEART AND SOUL...



**STERNOTOMY
DOUBLE WIRES
STAINLESS STEEL**

A perfect tool results in a perfect outcome

Innovative Sternotomy Wires from FUMEDICA



Advantages

- high degree of break resistance at twisting point
- large contact area at sternum
- reverse cutting needle ($\frac{3}{8}$ circle) for optimal puncture
- ease of use due to twister

Double Wires made of stainless steel, needle DS 48, sterile packed

wire size USP	wire length cm	wires per pouch	no. pouches per box	order number
7	2 x 30	1	12	DS4-770-1
7	2 x 30	2	12	DS4-770-2
7	2 x 30	3	12	DS4-770-3

Double Wires made of titanium, needle DS 48, sterile packed

wire size USP	wire length cm	wires per pouch	no. pouches per box	order number
7	2 x 30	1	12	TDS4-770-1
7	2 x 30	2	12	TDS4-770-2
7	2 x 30	3	12	TDS4-770-3

Twister device for twisting the double wire



order number

S4-100

Germany / International

Fumedica Medizintechnik GmbH

Lotzenäcker 17

D-72379 Hechingen

Phone +49 (0)7471 920 00 50

Fax +49 (0)7471 920 00 52

Internet www.fumedica.de

E-Mail fumedica@fumedica.de

Switzerland

Fumedica AG

Luzernerstrasse 91

CH-5630 Muri

Phone +41 (0)56 675 91 00

Fax +41 (0)56 675 91 09

Internet www.fumedica.ch

E-Mail fumedica@fumedica.ch

FUMEDICA 